

SAFETY 101 | Disinfect your devices

Five steps to disinfection¹

The digital devices you handle all day need to be disinfected just like other hospital equipment to keep you, your patients, and your team members safe. Here's how.



Use a non-abrasive, lint-free cloth for all steps. No paper towels!



1 Shut down the device and disconnect power or use the Lenovo Quick Clean application.



2 Wash hands for at least 20 seconds and dry thoroughly.



3 Apply an approved cleaning agent to the cloth (not onto the device) and completely wring excess liquid from the cloth.



4 Wipe all surfaces and allow the device to dry.



5 Reconnect power and turn on the device.

LENOVO QUICK CLEAN APPLICATION²

Lenovo products, like the ThinkPad T14 Gen 3 powered by Intel vPro[®] with 12th Gen Intel[®] Core™ processor, can be equipped with the Lenovo Quick Clean application. Lenovo Quick Clean will lock the touchscreen, keyboard, and mouse display automatically to allow for proper sanitation of the device between patients without shutting down the system. Customize your timer or use the 2-minute preset to lock out input.

DISINFECTING AGENTS

SAFE TO USE

- Hydrogen peroxide wipes (such as Virox Accel wipes w/ a maximum of 0.5% hydrogen peroxide)
- Alcohol-based wipes (such as PDI Germicidal or B.Braun Meliseptol Wipes sensitive)



NEVER USE chlorine or bleach

APPROXIMATELY

98% & **84.5%**
of healthcare workers own a mobile phone bring them to work every day.³

A bacterial and pathogen contamination rate of

94.3%
was observed on mobile phones and hands of healthcare workers.⁴



DON'T FORGET TO CLEAN ACCESSORIES!

For non-Lenovo devices, check with the original manufacturer's recommendations for cleaning and disinfecting of external keyboards, mice, cords, earbuds, and handsets.

HOW OFTEN?


Ideally, you should disinfect mobile devices before and after every patient encounter that involves contact — as well as before and after your shift.



PREVENTION PRACTICES Can reduce healthcare-associated infections by as much as 70%.⁵

LENOVO THINKPAD T14[®] GEN 3

This healthcare-ready workhorse features improved audio and new video cameras for virtual health conferencing, single sign-on enabled with RFID/NFC integrated badge reader, and heightened physical privacy with an electronic webcam shutter and optional privacy glass.



For more information on cleaning your PC, go to www.Lenovo.com/cleanpc

For more information on Lenovo Healthcare offerings, go to www.Lenovo.com/Health

Sources:

¹ <https://support.lenovo.com/us/en/solutions/ht035676>

² <https://support.lenovo.com/us/en/downloads/ds540666-lenovo-quick-clean-for-windows-10-thinkpad>

³ <https://www.sciencedirect.com/science/article/pii/S2590088921000147>

⁴ <https://www.sciencedirect.com/science/article/pii/S2666991921001160>

⁵ <https://aricjournal.biomedcentral.com/articles/10.1186/s13756-022-01125-8>



Smarter
technology
for all

Lenovo